

# ProBlot™ Rocker 25 and 25XL

Axon Labor Technik GmbH  
Im Braumenstück 16 - D-67659 Kaiserslautern  
Verkaufsbüro Tel. 06223/95 41 190 Fax: 06223/95 40 714  
info@axon-labortechnik.de

The Labnet compact rockers, ProBlot 25 and 25XL, are ideal for staining, hybridization, and general mixing. Their basic design and economical price makes them a great value for any laboratory.

Rocking speed is adjustable to provide the gentle motion required for staining fragile gels as well as the more vigorous action necessary for washing blots. The rocking movement is fixed at 7°.

The platform surfaces are covered with non-slip rubber pads to prevent sample containers from moving during operation. Elastic tie downs are available to secure bags, boxes, tubes and other containers in place. The ProBlot 25 platform measures 26 x 20 cm, while the ProBlot 25XL platform is significantly larger at 30 x 30 cm. For increasing the work area without sacrificing valuable bench space, double platform models are available.

## Features

- ▶ Can be used safely in temperature controlled environments up to 65°C
- ▶ Low profile and compact size make them perfect for use in incubators



S2025-D-B

## Specifications

Speed Range	5 - 120 tilts/min.
Timer	0 - 2 hr. or continuous
Motion	Rocking
Tilt Range	Fixed 7°
Clearance Between Platforms	3.5 in. (9 cm)
Maximum Load	11 lbs. (5 kg)
Ambient Operating Range	+4°C to 65°C
Dimensions (W x D x H)	
ProBlot 25	11.4 x 12.6 x 5.7 in. (29 x 32 x 14.5 cm)
ProBlot 25XL	12.8 x 14.2 x 5.7 in. (32.5 x 36 x 14.5 cm)
Weight	11.7 lbs. (5.3 kg)
Electrical	120V, 50/60Hz or 230V, 50/60Hz



S2025-XLD-B

## Ordering Information

Cat. No.	Description
S2025-B	ProBlot Rocker 25 with single platform (26 x 20 cm), 120V
S2025-B-230V	ProBlot Rocker 25 with single platform (26 x 20 cm), 230V
S2025-D-B	ProBlot Rocker 25D with double platform (26 x 20 cm), 120V
S2025-D-B-230V	ProBlot Rocker 25D with double platform (26 x 20 cm), 230V
S2025-XL-B	ProBlot Rocker 25XL with large single platform (30 x 30 cm), 120V
S2025-XL-B-230V	ProBlot Rocker 25XL with large single platform (30 x 30 cm), 230V
S2025-XLD-B	ProBlot Rocker 25XLD with large double platform (30 x 30 cm), 120V
S2025-XLD-B-230V	ProBlot Rocker 25XLD with large double platform (30 x 30 cm), 230V

## Accessories

Cat. No.	Description	Qty/Cs
S2025-T	Elastic tie downs	4