

LABNET INTERNATIONAL, INC.

MPS 1000tm Mini Plate Spinner

Finally, a compact and economical centrifuge for PCR plates

Labnet
Labnet International, Inc.



- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges

Patent Pending

SPECIFICATIONS

Capacity:	2 PCR Plates
Rotor Speed:	2500 rpm
G-Force:	500 xg
Rotor:	Fixed Vertically
Dimensions (WxDxH):	7.5 x 8.2 x 7.2 in. 19 x 21 x 18.3 cm
Weight:	6 lbs / 3 kg
Electrical:	120V or 230V, 50-60 Hz



ORDERING INFORMATION

Catalog No.	Description
C1000*	Mini Plate Spinner centrifuge with two position rotor

* To order this product in 230V, please add (-230V) to the item number.

Ihr Fachhändler:
www.axon-labortechnik.de

Axon
Labortechnik

Im Braumenstück 16 - D-67659 Kaiserslautern
Verkaufsbüro: 06222/93 56 80 Fax: 06222/93 56 66
eMail: info@axon-labortechnik.de